

AcrySof® IQ Toric takes precise astigmatism correction to a whole new level.





Precisely where you need to be.

Precise Astigmatism Correction

The AcrySof[®] IQ Toric IOL reduces astigmatism, for increased spectacle-independent distance vision and high patient satisfaction^{1,2}.

Your cataract surgery patients with astigmatism depend on you to restore their vision, and AcrySof® IQ Toric is the IOL you can trust to confidently provide precise astigmatism correction.

REDUCTION OF RESIDUAL REFRACTIVE CYLINDER



63% of patients implanted achieved \leq 0.50 diopters of residual refractive cylinder. 87% achieved \leq 1.00 diopters¹.

†† AcrySof® Single-Piece (SA60AT)

IMPROVED UNCORRECTED DISTANCE VISUAL ACUITY

225°

180°





Unparalleled Rotational Stability

The AcrySof[®] Single-Piece platform makes the difference.

Proven biomechanics and biomaterial ensure minimal rotation - less than 4° average rotation 6 months after implantation^{1,2}.

- > STABLEFORCE® haptics keep AcrySof® IQ Toric highly stable and centered in the capsular bag²
- > Flexible haptic design provides optimal placement in capsular bag, regardless of size²
- > AcrySof[®] lens material binds to fibronectin, ensuring adhesion to the anterior/posterior capsule³

IMPACT OF ROTATION **ON CORRECTION**



Generally, for every degree of IOL rotation, 3.3% of lens cylinder power is lost. A complete loss of cylinder power can occur with a rotation of 30° or greater².

ENS AXIS ORIENTATION ative vs 6 Months Postoperative)



81.1% of patients were $\leq 5^{\circ}$ of intended axis², and 971% were < 10° of intended axis¹



AcrySof® Aspheric Technology

Reduced Spherical Aberration

AcrySof® IQ Toric is designed with negative spherical aberration** to compensate for the positive aberration of the average cornea, which reduces both spherical and total higher order aberrations, for enhanced visual performance⁴.

-

SPHERICAL AND TOTAL HIGHER ORDER ABERRATIONS 90–120 DAYS AFTER 2ND EYE IMPLANT



A statistically significant reduction in both spherical and total higher order aberrations⁴.

- ★ Differences favor AcrySof® IQ IOL overall and at each visit (p<0.0001).</p>
- tt AcrySof® Single-Piece (SA60AT)

Increased Contrast Sensitivity

Engineered to improve contrast sensitivity in low-light conditions⁴, the aspheric design of AcrySof® IQ Toric plays a vital role in image quality.**

CONTRAST SENSITIVITY* IN MESOPIC CONDITIONS 90–120 DAYS AFTER 2ND EYE IMPLANT



AcrySof® IQ IOL showed statistically significant improvement⁴ in mesopic contrast sensitivity over the control lens in situations with and without glare at 6 cycles per degree (cpd) using the Vector Vision CSV-1000[§].**

- * Contrast sensitivity was measured using Vector Vision CSV-1000.
- § At 3 cpd, there was no significant change in mesopic contrast sensitivity.

MODULATION TRANSFER FUNCTION (MTF) VALUES 5 mm APERTURE⁵



AcrySof® IQ Toric IOL has similar performance to the AcrySof® IQ IOL in mesopic conditions.

** The effects of this aspheric design have been clinically assessed on AcrySof® IQ IOL Model SN60WF. See Directions for Use.

Delivers Exceptional Image Quality

Improved Functional Vision

HI I

Functional vision is an important consideration for your patients with astigmatism. When it comes to object detection and identification, **seconds can matter** — particularly in low visibility conditions, like night driving.

Offering clinically and statistically significant improvement in night driving test conditions⁴, the aspheric design of the AcrySof® IQ IOL provides **improved** functional vision in challenging, low-visibility environments.**

AcrySof® IQ FUNCTIONAL VISION STUDY4



AcrySof® IQ patients had an average increase of 130+ feet (versus the control lens) in which to stop after identifying a warning sign⁴.

Improved functional vision in challenging environments:

- > Patients were tested in night driving conditions that simulated city and rural settings under normal, glare and fog conditions.
- > Patients were asked to detect and identify various types of visual targets.

AcrySof® IQ IOL performed functionally better in 34 of 36 conditions tested:

- > Performance in 12 of these conditions was statistically significant, with the greatest advantage observed in detection and identification of city pedestrians and rural warning signs under glare and fog conditions.
- > AcrySof® IQ IOL improved functional vision (versus the control lens) by giving patients greater time to take appropriate action.



Expanded Range of Cataract Patients

With a wide range of cylinder powers, AcrySof® IQ Toric is designed to accommodate a variety of cataract patients with astigmatism, including those with lower levels of astigmatism.

ESTIMATED DISTRIBUTION OF PREOPERATIVE CYLINDER

90°

270°

135°

225°



*Based on average pseudophakic human eye.

0.75 to

1.50 D

1.50 to

2.00 D

2.00 D and up

Recommended

Corneal Astigmatism

Correction Range

www.AcrySofToricCalculator.com

AcrySof® IQ Toric Calculator

To ensure the best outcomes, you need a sound surgical plan. And precise surgical planning has never been easier, thanks to the AcrySof® IQ Toric Calculator. Easy input, powerful output—the AcrySof® IQ Toric Calculator is an innovative tool to help you perform your best.

Easy Input

- > Patient data
- > Keratometry
- > IOL spherical power
- > Surgically induced astigmatism
- > Incision location

Plase enter the pre-op information for the patient. Burgeon Name	m		and the second se		
Control Place enter the pre-op information for the patient. Surgeon Name	ACROCK	2			
Plane enter the pro-op information for the patient. Plane enter the pro-op information for the patient. Plane enter the pro-op information for the patient. Plane takes Pla	TORIC	Charles of Ch			
Storgen -0* Patient Name -0* Additional Patient Information Patient Information -0% K violatin 0 Dopter C stillenster File Selection 0 Dopter C stillenster File X ** 130* # File X asis e* - 118* 100* Steep Axis e* - 118* Dis Spherical Perser (fr-10) -0 +.8 D - 34.8 D Steep Gilly Include (SiA) -118* Steep Axis Dis Spherical Discult -2.880 Discision Lection (L) -2.880 Continue -2.880	Current Marrie	Please	e enter the pre-op inform	nation for the patient.	
Patient Name Additional Additional	Surgeon Name			90*	
Additional Additional Partian (LD, Case, ed.) Exestion C op (mphy) C os (Left) K Notation C Diopter C Hitmeter Plat A asis e*~18* Steep K 15.880 ~ 58.880 # Flat Asis e*~18* Steep K 15.880 ~ 58.880 # Steep Asis e*~18* Steep Asis e*~18* Dis Specifical - 2.880 Totsion - 2.880 Continue Continue Continue	Patient Name	1		and the second se	
Eye selection C op (mph) C op (mph) C op (mph) C op (mph) K Notation P Dopter C Hillmarker 180- 180- Flat. Asis e** - 118* 180- 180- 120- Steep Axis e** - 118* 120- 120- 120- Districtal - +.8 D - 34-8 D 120- 120- Steep Axis e* - 118* 118- 120- 120- Steep Axis e* - 2880 14- 120- 120- Incision - - 2.880 120- 120- Incision - - 120- 120- 120-	Additional Patient Information (I.D., Case, etc.)	-		135* 55*	
K Notation P Dupter 1 elimenter File K 15, 5400 - 58, 5400 # File Kais # - 188*1 Steep X 53, 580 - 58, 5800 Steep X 53, 580 - 58, 5800 Steep Axis # 18* Distribution 348.0 Swerge (-101) = 4, 80 - 34, 80 Swerge (-101) = 4, 80 - 34, 80 Swerge (-101) = 4, 80 - 34, 80 Distribution - 2, 580 Disclose - 2, 580 Lexation (IL) @ 383*	Eye Selection	C DD (Right) Cos (Left)	180%0*	
Flat K 15.480 - 58.480 # Flat Asis 4* - 118* Steep K 53.480 - 58.480 # Steep Asis - 218* IOL Spherical - 30 - 14.80 Surge Catly - 260 Discott - 2.860 Indicated - 2.860 Incision - 2.860 Lection (L) - 0* - 2880 Continue - 2.860	K Netation	C Diopter	C Hillimeter		
# Flat Axis #* ~ 18* 223* 315* Steep K 15.48D ~ 58.48D 200* 210* # Steep Axis #* ~ 18* 210* 210* # Steep Axis #* ~ 18* 11* 210* Steep Axis #* ~ 18* 11* 11* Steep Axis	Flat K	[35.00D ~ 50.00D		
Steep K 15.480 - 58.480 2/0* Ø Steep Axis 0* - 180* 2/0* Ø Steep Axis 0* - 180* 2/0* Steep Axis 0* - 180* 5/0* Steep Axis - 2.40 Flat Axis Steep Axis - 2.400 Flat Axis Incision - 2.400 Incision Iscation (II) - 2.400 Excation (II)	@ Flat Axis	[0*~180*	225* 315*	
ID: Spherical	Steep K	1	35.00D ~ 50.00D	270*	
101. Spherical	@ Steep Axis	[0*~180*		
Surgically Default (0.500) =.900 Included (S1A) - 2.860 Incision Lecation (1) Continue Continue	IOL Spherical Power (P-IOL)		6.8 D ~ 34.8 D	Steep Axis	
Location (L) Continue	Surgically Induced Astigmatism (SIA)		Default (0.500)	Incision	
Continue	Incision Location (IL)	1	0*~350*		
		Continu			

Powerful Output

- Recommended IOL model and spherical equivalent power
- > Optimal axis placement
- > Magnitude and axis of anticipated residual astigmatism





SPECIFICATIONS	
Model Number	SN6AT3, SN6AT4, SN6AT5
Optic Diameter	6.0 mm
Overall Length	13.0 mm
Optic Type	Biconvex Toric Aspheric Optic
IOL Powers (spherical equivalent diopters)	+6.0 to +30.0 D
IOL Cylinder Powers	1.50 D, 2.25 D and 3.00 D
Haptic Angulation	0 degrees (planar)
Haptic Configuration	STABLEFORCE® modified L haptic
Suggested A-Constant	119.0†
Refractive Index	1.55
Light Filtration	UV and Blue Light

†Provided as a guideline only.





References

- 1. Based on unilateral clinical study results (Models SA60T3, SA60T4, SA60T5). See package insert.
- 2. Data on file. Alcon, Inc.

Alcon

- Linnola RJ, Sund M, Ylönen R, Pihlajaniemi T. Adhesion of soluble fibronectin, laminin, and collagen type IV to intraocular lens materials. J Cataract Refract Surg. 1999;25:1486-1491.
- Results of a controlled, randomized, double-masked, multicenter, controlateral implant clinical study of the AcrySof® IQ IOL versus an AcrySof® Single-Piece IOL (SA60AT). See Directions for Use.
- Image quality was characterized by measuring MTF in a model eye that utilized a simulated cornea exhibiting typical adult human spherical aberration. Using the modified model eye, MTF measurements were made using both 3 and 5 mm apertures.

www.AcrySoflQToric.com



© 2009, Alcon, Inc. 3/09 TOR239

Precisely where you need to be.