icare HOME



IOP SELF-MEASURING
ANYWHERE, ANYTIME



THE 24HR TONOMETER

The easy-to-use Icare® tonometers revolutionize IOP measurements making it quick, effortless and effective. The Icare® HOME tonometer is designed for patients who need additional IOP monitoring according to their doctor's recommendation.

The device is based on a rebound measuring principle that eye care professionals have come to trust. A small probe touches the cornea and the impact duration, deceleration and other motion parameters are evaluated.

Icare® HOME does not require topical anesthetic or specialized skills for use. The HOME tonometer records IOP measurements and additional data which is not displayed to the patient. The healthcare professional downloads the measurement data via Icare LINK software.

KEY FEATURES



Icare® AMS
Automatic measuring sequence:
series and single mode
with one button



Icare® EasyPos
Intelligent positioning
assistant for the
correct alignment
of the tonometer



Icare® EyeSmart
Automatic OD/OS
recognition technology



- + REAL WORLD DATA
- + IDENTIFY IOP VARIATION AND SPIKES
- + VALUABLE ADJUNCT TO CLINIC

mmHg

CLINIC MEASUREMENTS (12 months)

- IOP monitoring in the clinic
- 24 hour IOP self-measuring with Icare HOME

mmHg







Icare HOME (24 HOURS)



ICARE USA - VISIONARY IN VISION

Icare USA is a subsidiary of Icare Finland, the original developer and global leader in handheld tonometry. Our patented technology (more than 20 patents and patent applications) combined with ISO 13485 certified quality system has made us a respected player in our field of expertise.

The advanced Icare® product line offers reliable, precise and reproducible accuracy when measuring intraocular pressure in both experienced and inexperienced hands. Icare tonometers available:

Icare[®] HOME, Icare[®] ic100, Icare[®] TA01i

Made in Finland.

ICARE HOME APP





ICARE USA

4208 Six Forks Road, Suite 1000, Raleigh, North Carolina, 27609 USA www.icare-usa.com

PROBE ORDERS: 888-422-7313 (Toll free)
SERVICE CENTER: 888-389-4022 (Toll free)
EMAIL: info@icare-usa.com